


LIFELONG[®]

SYRINGES



LIFELONG 1mL SYRINGE:
PRECISE & RELIABLE PERFORMANCE WITH
VARIOUS OPTIONS TO CHOOSE FROM

1 mL | SYRINGE
PREMIUM
MATRIX



LIFELONG 1 mL syringe for precise dosage preparation and injection. It offers wide range that provides solution for your different clinical needs. Ergonomic 3 parts design with minimal deadspace as per international standards which ensure accuracy in medication.



CATALOGUE

1 mL | SYRINGE











NEEDLE

- Made of high quality stainless steel
- Available in various lengths and diameters
- Ultra sharp long bevel, siliconized needle offers smooth penetration
- Siliconized needle surface further enhance the experience
- Color coded hub as per international standard for easy size identification

BARREL AND PLUNGER

- Transparent barrel with prominent graduation for easy reading
- Made of high quality Medical Grade Polypropylene
- Double seal gasket assembly to prevent leakage
- Wide finger grips allows easy application
- Secure plunger backstop to prevent accidental plunger withdrawal
- Specially engineered gasket spike ensures Low Dead Space

NEEDLE RANGE

S.No.	Needle Size		International Standard of Colour code
1	29G X 1/2"	0.33MM X 13MM	
2	27G X 1/2"	0.40MM X 13MM	
3	26G X 1/2"	0.45MM X 13MM	
4	25G X 1/2"	0.50MM X 13MM	
5	25G X 5/8"	0.50MM X 16MM	
6	25G X 3/4"	0.50MM X 19MM	
7	25G X 1"	0.50MM X 25MM	
8	25G X 1 1/2"	0.50MM X 38MM	
9	23G X 1"	0.60MM X 25MM	
10	23G X 1 1/2"	0.60MM X 38MM	

PACKAGING

S.No.	Product Desc.	Pieces/Box	Boxes/Case	Total Pcs/Case
1	1mL Syringe	100	30	3000



PVC
FREE

LATEX
FREE

DEHP
FREE

Lifelong Meditech Pvt. Ltd.

Factory & Corporate Office:

Plot No. 18, Sector 5, IMT Manesar,
Gurugram - 122050 (Haryana), India
Tel: +91-124-4406600 | Fax: +91-124-4406699

Follow us on:     

www.lifelongmed.com

Registered Office:

5th Floor, Unit No. 508, DLF South Court,
Saket District Centre, New Delhi - 110017, India
Tel: +91-11-26017063-64, 46017700 | Fax: +91-11-46017799

Email: info@lifelongindia.com

Customer Care Email: lmcustomer@lifelongindia.com

CIN No.: U85110DL2002PLC116251