

LIFELONG[®]

DISPOSABLE SYRINGES & NEEDLES



CATALOGUE

ABOUT THE COMPANY

Lifelong Group Overview

- ✓ More than **USD 175 Million Lifelong group** with interests ranging from auto-component and medical devices manufacturing to education.
- ✓ Group's flagship auto-components business is an award winning supplier to international OEMs such as GM and Hero Motocorp.
- ✓ Group has **3 plants in north India** located at Manesar, Dharuhera and Haridwar.



About Lifelong Meditech

- ✓ Started in 2005 to provide world class medical devices.
- ✓ In 17 years it has become one of the largest and fastest growing medical devices company in India.
- ✓ Supplier to prestigious government and private medical institutions across India and key International markets.
- ✓ Within a short span of time it has become the largest exporter of disposable syringes from India and commands 50% shares of indian disposable syringe exports.
- ✓ It has annual **installed capacity of 1.6 billion syringes** & produces disposable syringes ranging from 1 ml to 60ml.
- ✓ It is the first Indian company to bring world class Patented third generation technology fully **safe syringes** to India.

We welcome our customer to visit our factory.
Please call +91-124-4406600 to help us organize the same.


State of the art Manufacturing facility



10,000 Class Printing, Assembly & packing Area



High Quality Injection Molding Machine



Fully Automated Centrally Controlled Sterilization Facility


WORLD CLASS QUALITY



QUALITY POLICY

We are committed to:


- Providing quality products at optimum cost.
- Implementing continual improvement.
- Maintaining an effective quality management system which complies with regulatory requirement.
- Ensuring customer satisfaction and meeting their requirements on sustained basis.



QUALITY OBJECTIVES:

- Achieve customer satisfaction through timely & cost effective supply of products of consistent quality.
- Achieve continual improvement through:
 - Employee training and motivation
 - Productivity improvement
- Comply with all applicable Indian and International regulations.


All Product Batches are tested under statistical process in the World Class Laboratory of Company.



STANDARD SYRINGES: LUER SLIP & LUER LOCK

LIFELONG manufactures a full range of hypodermic syringes to match a variety of clinical needs. Our state-of-the-art manufacturing assures products of consistently high quality that also meet the most rigorous clinical requirements.

Syringes are made of non-toxic, medical grade virgin polypropylene and are latex free and pyrogen free.



- The tip of the barrel has a 6% lower taper as per ISO: 594
- Exceptional clarity for precise volume reading.
- High contrast black graduation with permanent marking ensures clear readability.
- Safe plunger made of Milky white polypropylene and has suitable protection to avoid accidental spill over.
- Gasket is made of highly inert TPE – thermoplastic elastomer in line with global practice and is compatible with any medication.
- Needle is made of cannula imported from world class manufacturer and has ultra thin wall permitting rapid and efficient drug delivery.
- Automated manufacturing process to ensure uniform quality.
- Sterilization with ETO gas using automated process as per ISO 11135.
- In blister pack syringes, Swedish medical grade paper and a good LD is used for better sealing and peel off.
- The design of the product is fully optimised in terms of cost and it complies with all ISO requirements.



NEW

Catheter Tip Syringe - 50ml



STANDARD SYRINGES

Manufactures of full range of hypodermic syringes to match a variety of clinical needs. Our state -of-art manufacturing assures products of consistently high quality that also meet the most rigorous clinical requirements.

PRODUCT FEATURES & QUALITY:

- Medical grade polypropylene
- Automated manufacturing process
- Consistent high quality, latex free Syringes
- Transparent barrel for prominent Reading

Available Sizes : 1, 2, 3, 5, 10, 20 & 50 ml



LIFELONG®

MATRIX

With & without Needle

LIFELONG®

PREMIUM

With Mounted Needle




HYPODERMIC NEEDLES

Lifelong hypodermic needle with top specification and raw material to ensure highest level of safety and efficiency of use for medical staff, minimal discomfort for patients. Needles are manufactured in compliance with ISO 7864 and it features the ultra sharp stainless steel through its optimized automated tipping technology for minimal penetration force. This contributes to preventing tissue injury and more commonly patient discomfort. Additionally LIFELONG 100% automated online inspection using cameras and sensors ensures total product quality.

PRODUCT FEATURES & QUALITY:

- Manufactured as per ISO 7864.
- Imported needle from world class manufacturer.
- Hub colour code for accurate, quick identification as per ISO 6009.
- Hub translucency provides confirmation of flashback.
- Needle wall is ultra thin permitting rapid and efficient drug delivery.
- Exceptionally sharp and three facet lancet point able to penetrate tissue smoothly and with minimal pain.
- Needles have ultra thin wall permitting rapid and efficient drug delivery.



Product Range

Gauge	Hub Color		Needle Length per gauge
18G (0.9mm)		Pink	1", 1 1/2"
19G (1.1mm)		Cream	1 1/2"
20G (0.9mm)		Yellow	1 1/2"
21G (0.8mm)		Green	1", 1 1/4", 1 1/2"
22G (0.7mm)		Black	1", 1 1/4", 1 1/2"
23G (0.6mm)		Blue	1", 1 1/4", 1 1/2"
24G (0.55mm)		Purple	1", 1 1/4"
25G (0.5mm)		Orange	5/8", 1", 1 1/2", 3/4"
26G (0.45mm)		Brown	1/2", 1", 3/8", 1 1/2"


Product Available on Request

Gauge	Hub Color		Needle Length per gauge
27G (0.4mm)		Grey	1/2"
28G (0.36mm)		Green	1/2"
29G (0.33mm)		Red	1/2"
30G (0.3mm)		Yellow	1/2"



CERTIFICATES & CREDENTIALS

- 🌐 Largest Syringe Exporter from India
- ★ Star Exporter Award from Indian Govt.
- 🌐 Exports to more than 40 Countries
- 🩺 Annual installed capacity of 1.6 billion syringes



LIFELONG[®]

Lifelong Meditech Pvt. Ltd.

Factory & Corporate Office:
Plot No. 18, Sector 5, IMT Manesar,
Gurugram - 122050 (Haryana), India
Tel: +91-124-4406600 | Fax: +91-124-4406699

Follow us on:     

www.lifelongmed.com

Registered Office:
5th Floor, Unit No. 508, DLF South Court,
Saket District Centre, New Delhi - 110017, India
Tel: +91-11-26017063-64, 46017700 | Fax: +91-11-46017799

Email: info@lifelongindia.com
Customer Care Email: imlcustomercare@lifelongindia.com
CIN No.: U85110DL2002PLC116251