



# LIFELONG<sup>®</sup>

## HYPODERMIC SAFETY NEEDLE

### SafeT<sup>®</sup>needle

LIFELONG takes pride in offering a new Safety product - Safety Needle. This needle helps in avoiding accidental needle pricks and infections risk. This new product will complement our existing range of Safety syringes.

**SafeT** is a sterile, disposable, hypodermic needle for injection administration with inneversible safety mechanism. It is easy to use. It is non-toxic, non-pyrogenic & sterilized by Ethylene Oxide. The individual blister packing, with needle hub on the side of peel open make easy assembling it on the syringe. Luer Slip & Luer Lock suitable for all kind of standard syringes.



CATALOGUE

# Features and benefits:










## Safety Needle Features:

- Protects from accidental pricks
- Provides irreversible, easy activation mechanism
- Suitable for any conventional Luer

## Technical Features:

- Stainless steel needle tube with triple bevel, diameter & length as per ISO-9626
- Polypropylene protective cap
- Safety system for needle cover in PP in the same color as the cone, according to International standards ISO 6009



AVAILABLE NEEDLE SIZES								
								
Pink	Beige	Yellow	Green	black	Blue	Purple	Orange	Brown
18G	19G	20G	21G	22G	23G	24G	25G	26G

## Lifelong Meditech Pvt. Ltd.

### Factory & Corporate Office:

Plot No. 18, Sector 5, IMT Manesar,  
Gurugram - 122050 (Haryana), India  
Tel: +91-124-4406600 | Fax: +91-124-4406699

Follow us on:     

[www.lifelongmed.com](http://www.lifelongmed.com)

### Registered Office:

5th Floor, Unit No. 508, DLF South Court,  
Saket District Centre, New Delhi - 110017, India  
Tel: +91-11-26017063-64, 46017700 | Fax: +91-11-46017799

Email: [info@lifelongindia.com](mailto:info@lifelongindia.com)  
Customer Care Email: [imlcustomercare@lifelongindia.com](mailto:imlcustomercare@lifelongindia.com)  
CIN No.: U85110DL2002PLC116251